

sf aid

ΙΑΤΡΙΚΟΣ ΕΞΟΠΛΙΣΜΟΣ
ΑΝΑΛΩΣΙΜΑ

ΑΠΟΣΤΕΙΡΩΣΗ



ΥΓΡΟΣ ΚΛΙΒΑΝΟΣ HYDRA EVO AUTOCLAVE - 15 It



Τιμή: Επικοινωνήστε μαζί μας

Technical Specifications

- Stainless steel chamber: \varnothing 24.5 x 32 cm
- External size: 44 x 38 x 50 cm (W x H x D)
- Net weight: 35 kg
- 2 sterilization cycles: 121°C and 134°C
- Operating Voltage: 230V \pm 10% - 50/60Hz
- Power consumption: 1700 W
- Class N
- Made in Italy - Warranty: 1 year

- Cycle: Not wrapped 121° C
- Materials:
 - Fragile solids
 - Rubber solids
 - Glass
 - Turbines (if indicated by producer)
- Sterilization time: 20 min
- Drying time: 5 min
- Max load: 2.8 kg (700 g load in any tray)
- Working pressure: 1.2 - 1.3 bars

- Cycle: Not wrapped 134° C
- Materials:
 - Stainless metal solids
 - Turbines (if indicated by producer)
- Sterilization time: 10 min
- Drying time: 5 min
- Max load: 2.8 kg (700 g load in any tray)
- Working pressure: 2.1 - 2.2 bar

ΥΓΡΟΣ ΚΛΙΒΑΝΟΣ HYDRA EVO AUTOCLAVE - 18 It



Τιμή: Επικοινωνήστε μαζί μας

Technical Specifications

- Stainless steel chamber: \varnothing 24.5 x 40 cm (Plus)
- External size: 44 x 38 x 50 cm (W x H x D)
- Net weight: 35 kg
- 2 sterilization cycles: 121°C and 134°C
- Operating Voltage: 230V \pm 10% - 50/60Hz
- Power consumption: 1700 W
- Class N
- Made in Italy - Warranty: 1 year

- Cycle: Not wrapped 121° C
- Materials:
 - Fragile solids
 - Rubber solids
 - Glass
 - Turbines (if indicated by producer)
- Sterilization time: 20 min
- Drying time: 5 min
- Max load: 2.8 kg (700 g load in any tray)
- Working pressure: 1.2 - 1.3 bars

- Cycle: Not wrapped 134° C
- Materials:
 - Stainless metal solids
 - Turbines (if indicated by producer)
- Sterilization time: 10 min
- Drying time: 5 min
- Max load: 2.8 kg (700 g load in any tray)
- Working pressure: 2.1 - 2.2 bar

ΞΗΡΟΣ ΚΛΙΒΑΝΟΣ ΤΑΥ STERIL AUTO DRY HEAT 18lt



Technical Specifications

- Inside dimensions: 375 x 230 x h 200 mm
- Outside dimensions: 440 x 310 x h 380 mm
- Weight: 14 kg
- Operating voltage: 220V (110 V available on request)
- Power consumption: 600 W
- Fuses: 3.15 A
- Thermostat: Automatic capillary bulb type with self registration from 90 to 200 °C
- Timer: from 0 to 120 minutes
- Thermometer: high precision with capillary bulb
- Made in Italy

Τιμή: Επικοινωνήστε μαζί μας

Andromeda Vacuum XR Αυτόκλειστο -15 λίτρα - 230V



Technical Specifications

- External size: 45 x 39 x 54 cm (W x H x D)
- Chamber dimension: Ø 24.5 x 32 cm depth
- Weight: 50 kg
- Operating Voltage: 230V - 50/60 Hz
- Power consumption: 2000 W
- Made in Italy
- Class S
- Cycles:

Andromeda Vacuum is provided with 5 cycles of sterilization: 121° C and 134° C for unwrapped and wrapped instruments, and a "Flash" cycle for the fast sterilizations.

- Program 1 (Unwrapped instruments)
- Temperature: 134° C
- Max load: 2.8 Kg
- Pressure: 2.2
- Total time*: 33 min

- Program 2 (Unwrapped instruments)
- Temperature: 134° C
- Max load: 2.8 Kg
- Pressure: 2.2
- Total time*: 38 min

- Program 3 (Unwrapped instruments)
- Temperature: 121° C
- Max load: 2.8 Kg
- Pressure: 1.2
- Total time*: 42

- Program 4 (Unwrapped instruments)
- Temperature: 121° C
- Max load: 2.8 Kg
- Pressure: 1.2
- Total time*: 46

- Program 5 (Flash)
- Temperature: 134° C
- Max load: 2.8 Kg
- Pressure: 1.2
- Total time*: 16

* pre-heating - vacuum - cycle - drying in vacuum

Τιμή: Επικοινωνήστε μαζί μας

NUBYRA Αυτόκλειστο - 6 λίτρα - 230V



Τιμή: Επικοινωνήστε μαζί μας

Technical Specifications

- Anodized aluminium chamber: \varnothing 17x27 cm - 6.1l
- External size: 44 x 45 x h 35 cm
- Trays size:
 - 1 large 15 x 24 x 2 cm
 - 2 small 11 x 24 x 2 cm
- 2 sterilization cycles: 121°C and 134°C
- Maximum load: 2 kg
- Operating voltage: 230V - 50/60 Hz
- Max absorbed power: 1,300 Watt
- Net weight: 27 kg Class N
- Made in Italy - Warranty: 1 year
- User manual: English, Italian, (French, Spanish on request).

- Cycle: Not wrapped 121° C
- Materials:
 - Fragile solids
 - Rubber solids
 - Glass
 - Turbines if indicated by producer
- Sterilization time: 18 min
- Drying time:
 - 5 min short dry
 - 15 min long dry
- Max Load: 2.1 kg (700 g load in any tray)
- Working pressure: 1.1-1.2 bars

- Cycle: Not wrapped 134° C
- Materials:
 - Stainless metal solids
 - Turbines if indicated by producer
- Sterilization time: 6 min
- Drying time:
 - 5 min short dry
 - 15 min long dry
- Max Load: 2.1 kg (700 g load in any tray)
- Working pressure: 2.1-2.2 bars

ΞΗΡΟΣ αποστειρωτής "TAU 1.6 - 7lt"



Technical Specifications

- Inside dimensions: 245 x 115 x h 55 mm
- Outside dimensions: 310 x 180 x h 180 mm
- Weight: 5 kg
- Operating voltage: 220V (110 V available on request)
- Power consumption: 250 W
- Thermostat: Automatic capillary bulb type with self registration from 90 to 200 °C
- Timer: from 0 to 120 minutes
- Thermometer: high precision with capillary bulb
- Made in Italy.

Τιμή: Επικοινωνήστε μαζί μας



Αισχύλου 83, Νίκαια
Τηλ.: 211 4112707 - 6939 080534, 6931 981515
E-mail: info@sfaid.gr - Site: www.sfaid.gr